EXHIBIT A



Quotation

Date	Quote #
5/1/2015	20124514

Mailing Address: 1285 Wallis Street Eugene, OR 97402 Sales Representative: Erik Richter 916.747.0816 erichter@commfitness.com

Bill To:

Mendocino Coast Recreation and Park Distr C.V. Starr Community Center Ethan Newton 300 S Lincoln St Fort Bragg, CA 95437 Ship To:

C.V. Starr Community Center Ethan Newton 300 S Lincoln St Fort Bragg, CA 95437

Terms: Net 30 (OAC)

Qty Description Cost Total Item SciFit Step One Total Body Recumbent Stepper -SONE03 4,466.00 8,932.00 adjustable/swivel seat F/I Freight and inside Delivery shipped direct 695.00 695.00 TAXCA Sales Tax 8.625% 770.39 770.39

Quote valid for 30 Days. 50% Deposit and balance after Delivery. Net 30 Terms available upon Approved Credit

Total

\$10,397.39

Signature

Phone #	Fax#	E-mail
541-343-3790	541-343-3795	cmay@commfitness.com

Commercial Fitness Equipment Standard Terms and Conditions Apply





StepOne™ Product Specifications

Order Code: SONEØ1

FEATURES	LABRES SELECTION OF THE LOCAL PROPERTY OF TH
Smooth, linear movement with start and stop motion	Yes
Low starting resistance	6 watts starting resistance
Resistance range	20 levels adjustable in .1 increments, providing 200 levels of resistance
Neutral starting position	Pedals always come back to neutral starting position for easier entry and start up. Allows for entry from both sides.
Arm to leg movement ratio	1:1 Natural arm to leg movement ratio. Patterning similar to walking or climbing stairs.
User-defined range of motion	2" to 11.5" step length.
Comfortable, ergonomic body position	Yes
Total body workout	Yes
Large flywheel	Increased kinetic energy and smooth stepping motion.
Small footprint	5" shorter than other leading recumbent steppers.
Easy to transport	Integrated transport handle and wheels.
Heart rate monitoring	Compatible with most wireless technology.
Quiet belt drive	Yes
ACCESSIBILITY &	
Direct wheelchair access	Seat slides off easily for direct wheelchair access. ADA compliant.
Easy patient transfer	Adjustable swivel seat sits on a flat horizontal monorail and provides easy patient transfer.
Step-through seating	True adjustable step-through design, allows all users full access and ease of entry.
Step-through height	2"
Grab bar	Integrated for easy entry and support.
SEAT	DANS DESCRIPTION OF TRANSPORT
Swivel seat	Seat moves fore and aft and swivels 360 degrees, locking at each 90 degree angle.
Horizontal seat monorail	Seat adjusts forward and back on a flat monorail for easy and safe adjustment.
Removable seat	Yes, seat removes easily without hardware for direct wheelchair access
Horizontal seat adjustment range	18"
Optional bariatric seat	Yes - 600 lb. weight capacity. (Order code SONEØ2)
ARMS & GRIPS	LANGE DATE DE LES TRANSPORTES DE LA COMPANION
Arm length adjustment	: Adjust length of arms within a 12" range.
60 degree hand grip rotation	Position handles anywhere within the 60 degrees of hand grip rotation
Non-absorbent, chemical resistant hand grips	Easy to clean, durable, made from 100% recyclable material.
Easy adjustment levers	Yes

		-	 -	-		
-	- (0) -	n	Λ	6
			 _	_	$\overline{}$	-

Orthopedic footbeds with safety edge	Sturdy, oversized footbeds.
Heel cup	Yes
Foot straps	Yes

INTELLI-FIT™ CONSOLE

7" LCD full color touch screen display	Yes				
Easy operation	Touch screen with clear, simple screens and large, readable font.				
Audible beep	Yes				
Tactile buttons	Tactile buttons in	Tactile buttons in addition to touch screen.			
Resistance levels	200 levels of resi	stance - 20 levels adjus	stable in .1 increments.		
Display feedback	Time, total steps, steps per minute, stride length, calories, ascent, METs, watts, heart rate and level.				
Stride length	Stride length is accurately measured and displayed to help monitor range of motion and improvement.				
Ascent	Displays a true m	neasure of ascent (ft.) b	ased on stride length and steps.		
One touch quick start	Yes				
Cool down	Automatic cool down begins after each exercise is completed.				
Administrative mode	User defined set settings.	tings, metric or US unit	s of measure and other custom		
Download results to USB flash drive	Plug USB flash d	rive into the Fit-Key® sl	ot to save workout results in a		
CSAFE port and 8-volt DC power port	Yes				
USB charger	For phones, mp3 players and other electronic devices.				
Water bottle holder	Integrated				
Fan	Integrated dual fans with three speeds.				
Programs	Quick Start	Fit-Quik®	Constant Work (constant Watts)		
	Manual	Stress Test	Iso-Strength (constant RPM)		
	Heart Rate	7 Hill Profiles	Random (unlimited profiles)		

TECHNICAL SPECIFICATIONS

Maximum user weight capacity	450 lbs. standard seat; 600 lbs. bariatric seat.
User height range	Fits users 4'6" - 6"7".
Accuracy to EN 957-8 standard	Yes
Power Requirements	Self generating with auto recharge battery backup. Low voltage AC adapter included.
Unit weight	275 lbs.
Unit size	68"L x 30"W x 48"H
Manufacturing	Assembled in USA.

WARRANTY	TYPE ZAME DALLAM TULLANTAL OF
Parts	5 years
Labor	1 year
Frame	10 years
ENHANCEMENTS	
Fit-Key® software	Optional Fit-Key® software makes it easy to create specialized programs and track data.
Bariatric seat	Optional seat has 600 lb. user weight capacity, seat slides off easily for direct wheelchair access.
Entertainment options	Yes
ACCESSORIES	
Leg Stabilizer	Attaches at the calf where support is needed most. Natural pivot point at ankle. Durable steel construction.
Wheelchair Platform	Provides a sturdy foundation for wheelchair use.

Wrist support and enhanced grip - helps user maintain grip and

Transmits heart rate directly to console - works on all SCIFIT products.



SCIFIT Systems, Inc. reserves the right to change or alter specifications at any time.

©Copyright 2013, SCIFIT Systems, Inc. SCIFIT, Fit-Quik, and Fit-Key, are registered trademarks of SCIFIT Systems, Inc. Intelli-Fit, Heart Fit, Power Fit, StepOne, and Scientific Solutions for Fitness are trademarks of SCIFIT Systems, Inc. Fit-Key Technology and StepOne Recumbent Stepper are Patented.

StepOne Specs - 2013 - v2

Assist Gloves

Polar Chest Strap





StepOne™ with Bariatric Seat Product Specifications

Order Code: SONEØ2

EATURES	
Smooth, linear movement with start and stop motion	Yes
Low starting resistance	6 watts starting resistance
Resistance range Neutral starting position	20 levels adjustable in .1 increments, providing 200 levels of resistance Pedals always come back to neutral starting position for easier entry and start up. Allows for entry from both sides.
Arm to leg movement ratio	1:1 Natural arm to leg movement ratio. Patterning similar to walking or climbing stairs.
User-defined range of motion	2" to 11.5" step length.
Comfortable, ergonomic body position	Yes
Total body workout	Yes
Large flywheel	Increased kinetic energy and smooth stepping motion.
Small footprint	5" shorter than other leading recumbent steppers.
Easy to transport	Integrated transport handle and wheels.
Heart rate monitoring	Compatible with most wireless technology.
Quiet belt drive	Yes
CCESSIBILITY &	
CCE33IBIEITT G.	
Direct wheelchair access	Seat slides off easily for direct wheelchair access. ADA compliant.
Easy patient transfer	Seat adjusts forward and back on a flat horizontal monorail and provides easy patient transfer.
Step-through seating	True adjustable step-through design, allows all users full access and ease of entry.
Step-through height	2"
Grab bar	Integrated for easy entry and support.
EAT	DAMENTAL DATE TO A LOCAL TO A LOC
Bariatric seat	Oversized bariatric seat accommodates users up to 600 lbs.
Bariatric seat Horizontal seat monorail	Oversized bariatric seat accommodates users up to 600 lbs. Seat adjusts forward and back on a flat monorail for easy and safe adjustment.
•••••	Seat adjusts forward and back on a flat monorail for easy and safe
Horizontal seat monorail	Seat adjusts forward and back on a flat monorail for easy and safe adjustment.
Horizontal seat monorail Removable seat	Seat adjusts forward and back on a flat monorail for easy and safe adjustment. Seat removes easily without hardware for direct wheelchair access.
Horizontal seat monorail Removable seat Horizontal seat adjustment range	Seat adjusts forward and back on a flat monorail for easy and safe adjustment. Seat removes easily without hardware for direct wheelchair access.
Horizontal seat monorail Removable seat Horizontal seat adjustment range RMS & GRIPS	Seat adjusts forward and back on a flat monorail for easy and safe adjustment. Seat removes easily without hardware for direct wheelchair access. 18"
Horizontal seat monorail Removable seat Horizontal seat adjustment range RMS & GRIPS Arm length adjustment	Seat adjusts forward and back on a flat monorail for easy and safe adjustment. Seat removes easily without hardware for direct wheelchair access. 18" Adjust length of arms within a 12" range.

StepOne™ with Bariatric Seat Product Specifications (cont.)

	1 0
FOOT PEDA	

Orthopedic footbeds with safety edge	Sturdy, oversized footbeds.
Heel cup	Yes
Foot straps	Yes

Order Code: SONEØ2

INTELLI-FIT™ CONSOLE

7" LCD full color touch screen display	Yes				
Easy operation	Touch screen with clear, simple screens and large, readable font.				
Audible beep	Yes				
Tactile buttons	Tactile buttons	Tactile buttons in addition to touch screen.			
Resistance levels	200 levels of re	200 levels of resistance - 20 levels adjustable in .1 increments.			
Display feedback	•	Time, total steps, steps per minute, stride length, calories, ascent, METs, watts, heart rate and level.			
Stride length		accurately measured n and improvement.	and displayed to help monitor		
Ascent	Displays a true	measure of ascent (ft.) based on stride length and steps.		
One touch quick start	Yes				
Cool down	Automatic cool down begins after each exercise is completed.				
Administrative mode	User defined se settings.	ttings, metric or US u	nits of measure and other custom		
Download results to USB flash drive	Plug USB flash of the state of	drive into the Fit-Key [®]	slot to save workout results in a		
CSAFE port and 8-volt DC power port	Yes				
USB charger	For phones, mp3 players and other electronic devices.				
Water bottle holder	Integrated				
Fan	Integrated dual	fans with three spee	ds.		
Programs	Quick Start	Fit-Quik®	Constant Work (constant Watts)		
	Manual	Stress Test	Iso-Strength (constant RPM)		
	Heart Rate	7 Hill Profiles	Random (unlimited profiles)		

TECHNICAL SPECIFICATIONS

Maximum user weight capacity	600 lbs. (bariatric seat)
User height range	Fits users 4'6" - 6"7".
Accuracy to EN 957-8 standard	Yes
Power Requirements	Self generating with auto recharge battery backup. Low voltage AC adapter included.
Unit weight	275 lbs.
Unit size	68″L x 30″W x 48″H
Manufacturing	Assembled in USA.

StepOne™ with Bariatric Seat Product Specifications (cont.)

WARRANTY		
Parts	5 years	
Labor	: 1 year	
Frame	10 years	

ENHANCEMENTS

Fit-Key® software	Optional Fit-Key® software makes it easy to create specialized programs
	and track data.
Entertainment options	Yes

ACCESSORIES

Leg Stabilizer	Attaches at the calf where support is needed most. Natural pivot point at ankle. Durable steel construction.
Wheelchair Platform	Provides a sturdy foundation for wheelchair use.
Assist Gloves	Wrist support and enhanced grip - helps user maintain grip and position.
Polar Chest Strap	Transmits heart rate directly to console - works on all SCIFIT products.





Order Code: SONEØ2

Bariatric Seat 600 lb. user weight capacity

SCIFIT Systems, Inc. reserves the right to change or alter specifications at any time.

©Copyright 2013, SCIFIT Systems, Inc. SCIFIT, Fit-Quik, and Fit-Key, are registered trademarks of SCIFIT Systems, Inc. Intelli-Fit, StepOne, and Scientific Solutions for Fitness are trademarks of SCIFIT Systems, Inc. Fit-Key Technology and StepOne Recumbent Stepper are Patented.

StepOne Specs - SONE02 - 2013 - v1





StepOne™ with Height Adjustable Seat Product Specifications

Order Code: SONEØ3

FEATURES	
Smooth, linear movement with start and stop motion	Yes
Low starting resistance	6 watts starting resistance
Resistance range	20 levels adjustable in .1 increments, providing 200 levels of resistance
Neutral starting position	Pedals always come back to neutral starting position for easier entry and start up. Allows for entry from both sides.
Arm to leg movement ratio	1:1 Natural arm to leg movement ratio. Patterning similar to walking or climbing stairs.
User-defined range of motion	2" to 11.5" step length.
Comfortable, ergonomic body position	Yes
Total body workout	Yes
Large flywheel	Increased kinetic energy and smooth stepping motion.
Small footprint	5" shorter than other leading recumbent steppers.
Easy to transport	Integrated transport handle and wheels.
Heart rate monitoring	Compatible with most wireless technology.
Quiet belt drive	Yes
ACCESSIBILITY &	
Direct wheelchair access	Seat slides off easily for direct wheelchair access. ADA compliant.
Easy patient transfer	Seat adjusts up and down, forward and back, and swivels on a flat horizontal monorail and provides easy patient transfer.
Step-through seating	True adjustable step-through design, allows all users full access and ease of entry.
Step-through height	2"
Grab bar	Integrated for easy entry and support.
EAT	EASE DE DE DE DE DE DE DE DE
Height adjustable, swivel seat	Seat moves fore and aft, up and down, and swivels 360 degrees, lockin at each 90 degree angle.
Horizontal seat monorail	Seat adjusts forward and back on a flat monorail for easy and safe adjustment.
Removable seat	Seat removes easily without hardware for direct wheelchair access.
Horizontal seat adjustment reange	18"
Optional bariatric seat	Yes - 600 lb. weight capacity. (Order code SONEØ2)
RMS & GRIPS	
Arm length adjustment	Adjust length of arms within a 12" range.
60 degree hand grip rotation	Position handles anywhere within the 60 degrees of hand grip rotation
Non-absorbent, chemical resistant hand grips	Easy to clean, durable, made from 100% recyclable material.
Easy adjustment levers	Yes

OOT PEDALS	
Orthopedic footbeds with safety edge	Sturdy, oversized footbeds.
Heel cup	Yes
Foot straps	Yes
ITELLI-FIT™ CONSOLE	THE PARTY OF THE P
7"LCD full color touch screen display	Yes
Easy operation	Touch screen with clear, simple screens and large, readable font.
Audible beep	Yes
Tactile buttons	Tactile buttons in addition to touch screen.
Resistance levels	200 levels of resistance - 20 levels adjustable in .1 increments.
Display feedback	Time, total steps, steps per minute, stride length, calories, ascent, ME watts, heart rate and level.
Stride length	Stride length is accurately measured and displayed to help monitor range of motion and improvement.
Ascent	Displays a true measure of ascent (ft.) based on stride length and step
One touch quick start	Yes
Cool down	Automatic cool down begins after each exercise is completed.
Administrative mode	User defined settings, metric or US units of measure and other custor settings.
Download results to USB flash drive	Plug USB flash drive into the Fit-Key® slot to save workout results in a .txt file.

Yes

Integrated

Quick Start

Heart Rate

Manual

TECHNICAL SPECIFICATIONS

CSAFE port and 8-volt DC power port

USB charger

Fan

Programs

Water bottle holder

Maximum user weight capacity	450 lbs. (height adjustable seat)
User height range	Fits users 4'6" - 6"7".
Accuracy to EN 957-8 standard	Yes
Power Requirements	Self generating with auto recharge battery backup. Low voltage AC adapter included.
Unit weight	275 lbs.
Unit size	68"Lx30"Wx48"H
Manufacturing	Assembled in USA.

For phones, mp3 players and other electronic devices.

Constant Work (constant Watts)

Iso-Strength (constant RPM)

Random (unlimited profiles)

Integrated dual fans with three speeds.

Fit-Quik®

Stress Test

7 Hill Profiles

StepOne™ with Height Adjustable Seat Product Specifications (cont.) Order Code: SONEØ3

WARRANTY	
Parts	5 years
Labor	1 year
Frame	10 years
ENHANCEMENTS	
Fit-Key® software	Optional Fit-Key® software makes it easy to create specialized programs and track data.
Entertainment options	Yes
ACCESSORIES	TWATANDE MEDILANDE OF
Leg Stabilizer	Attaches at the calf where support is needed most. Natural pivot point at ankle. Durable steel construction.
Wheelchair Platform	Provides a sturdy foundation for wheelchair use.
Assist Gloves	Wrist support and enhanced grip - helps user maintain grip and position.
Polar Chest Strap	Transmits heart rate directly to console - works on all SCIFIT products.



SCIFIT Systems, Inc. reserves the right to change or alter specifications at any time.

©Copyright 2013, SCIFIT Systems, Inc. SCIFIT, Fit-Quik, and Fit-Key, are registered trademarks of SCIFIT Systems, Inc. Intelli-Fit, StepOne, and Scientific Solutions for Fitness are trademarks of SCIFIT Systems, Inc. Fit-Key Technology and StepOne Recumbent Stepper are Patented.

StepOne Specs - SONE03 - 2013 - v1